



GENiO™ User Training

GENiO™
powered by PDi

3 WAYS TO CLEAR PATIENT DATA WITH GENiO™

1

Terminals Tab

Select a terminal from the Terminals Tab, Click the "Clear User Data" Action icon to clear data cache. A message will show on your screen to notify you that all patient data has been erased.



2

On Device

On medTAB devices touch the "Clear Data" Icon at the bottom right of the screen. On A-Series devices, select the "Clear Data" icon or press the "0" button on the pillowspeaker or remote.



Clear Data

0-Clear Data

3

Power Off

When the device is powered down, all user data is automatically cleared.

