

DESIGN GUIDE

PDi TV SOLUTIONS TO DEPLOY EPIC MYCHART BEDSIDE TV

Provide a TV-centric experience to patients in the hospital

PDi
Communication Systems, Inc.



How To Qualify

- Talk to your Epic MyChart representative and have them reference SLG 7083467
- Bedside Subscriber
- Minimum Epic Version May 2024

PDi Android™-based Hospital TVs

- 75" - 65" - 55" - 43" Screen Sizes
- UL-Listed for Healthcare
- TV Mounts also available from PDi
- PDi TV Box for retrofit TVs

How to Deploy Bedside TV App

- PDi GENiO™ RDM cloud software subscription
- Network Access to whitelist GENiO servers

How Patients Control TVs

- Smart pillow speaker configured to your nurse call system
- Requires 4-way navigation, Back, Home, and OK/Enter buttons
- Recommended: Curbell Medical PDi-GENiO PS
- Submit front and back photos for compatibility

Link to 3rd Party Apps

- Ask PDi ProServices
- Analog, Digital, or IPTV Signal
- Virtual Care 2-way Video

Live Sites

- ✓ UT Health, TX
- ✓ Mercy Health, MD
- ✓ Altru Health, ND



Next Steps

1. Where to Buy – Contact PDi
2. Request Demo Equipment
3. Schedule GENiO Webinar
4. Configure Pillow Speaker
5. Pre-order TVs to Meet Deadlines



For more than 45 years, Ohio-based PDi has delivered premium reliability in patient-room TV technology. We bring Epic Bedside TV customers the all-inclusive technical knowledge needed to plan, deploy, and support successful smart room projects. From product design and infrastructure planning through installation and go-live, our team helps hospital staff manage every aspect of implementation with confidence.

Contact PDi
800.628.9870

Cyndi Kisenko ext 483
ckisenko@pdiam.com

Ric Torres ext 225
rtorres@pdiam.com

Rich Chuck ext 484
rchuck@pdiam.com

Jeff Conrath ext 223
jconrath@pdiam.com